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
Table 1 - Product and test parameters for extractables studies

	Sterile Connector	Capsule Filter
Liquid contact materials	Polycarbonate	Membrane: PVDF surface-modified with hydrophilic acrylate copolymer Support an drainage: Nonwoven polypropylene Capsule shell/filter hardware: Polypropylene
Contact surface area	Approximately 44 cm ²	Membrane: 1,500 cm ² Filter hardware: 500 cm ²
Extraction volume	1.5 L	1.5 L
Extracting fluid	Deionized water at 21°C; 18 megohm-cm resistivity; ethanol at 21°C	
Extraction system	Pumped recirculating liquid; PTFE pump and tubing	
Extraction time	Deionized water for 24 hours; ethanol for four hours	
Product pretreatment	Gamma irradiation at 50 kGy	

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Table 2 - Analytical methods used for assessing extractables

Analytical Method	Water Extract	Ethanol Extract
Total organic carbon (TOC) analysis	√	
pH analysis	√	
Conductivity analysis	√	
Ion chromatography (IC)	√	
Nonvolatile residue (NVR) analysis	√	√
Fourier-transform infrared (FTIR) spectroscopy	√	√
Ultraviolet (UV) spectrometry	√	√
Direct-injection gas chromatography with mass spectrometric detection (GC/MS)	√	√
Headspace GC/MS	√	√
Derivatization GC/MS	√	√
High-performance liquid chromatography with ultraviolet detection (HPLC/UV)	√	√
Liquid chromatography with mass spectrometric detection (LC/MS and LC/MS/MS)	√	√
Inductively coupled plasma mass spectrometry (ICP/MS)	√	√

Table 3 - Results of extractables studies on sterile connectors
(NA=not applicable; DL=detection limit)

Analytical Test Method	Solvent: Water		Solvent: Ethanol	
	Sterile connector	Negative Control	Sterile connector	Negative Control
Nonvolatile residue analysis	<0.05 mg	<0.2 mg (DL)	≤0.05 mg	<0.2 mg (DL)
Total organic carbon analysis	0.039 ppm	0.07 ppm	NA	NA
pH analysis (from a set of four connectors)	5.9	5.8	NA	NA
Conductivity analysis (from a set of four connectors)	1.9 μS/cm	1.0 μS/cm	NA	NA
Ion chromatography				
Acetate	<7.5 ppb ¹	<7.5 ppb ¹	NA	NA
Formate	<4.0 ppb ¹	<4.0 ppb ¹	NA	NA
GC/MS				
Volatiles (headspace)	No peaks	No peaks	No peaks	No peaks
Semivolatiles (direct injection)	No peaks ²	No peaks	No peaks	No peaks
Derivatization (BF ₃ /1-Butanol)	Benzoic acid: 0.039 ppm	Benzoic acid: 1.011 ppm	Benzoic acid: 0.047 ppm	Benzoic acid: <0.001 ppm
HPLC/UV	No peaks ²	No peaks	No peaks ³	No peaks
ICP/MS				
15 elements ⁴	<DL	<DL	<DL	<DL
Boron	0.397 ppb	0.08 ppb	<DL	<DL
Calcium	<DL	8.938 ppb	1.028 ppb	4.316 ppb
Sodium	0.162 ppb	0.07 ppb	<DL	47.64 ppb
Zinc	0.050 ppb	0.01 ppb	<DL	<DL

¹ Limit of detection: 7.5 ppb for acetate and 4.0 ppb for formate
² Extract also concentrated 50-fold (no peaks)
³ Extract also concentrated 10-fold (no peaks)
⁴ Li, Mg, Al, K, Ti, Cr, Mn, Fe, Co, Ni, Cu, Ag, Sn, Ba, and Pb

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Table 4 - Nonvolatile residues (NVR) from capsule filters

Samples	NVR
Water extract	4.1–6.8 mg ¹
Ethanol extract	50.8–60.9 mg ²
¹ After subtraction of negative control (0.3 mg)	
² After subtraction of negative control (0.6 mg)	

Table 5 - Identity and quantities of some compounds in ethanol filter extracts

Test	Compounds Identified	Concentration	
GC/MS on ethanol extract	2-Ethylhexanoic acid	0.56 ppm	
	1,3-Di-tert-butyl-benzene	0.52 ppm	
	2,4-Di-tert-butyl-phenol	0.12 ppm	
	Lauryl acetate	0.13 ppm	
	Lauryl acrylate	0.64 ppm	
Test	Fatty Acid Derivatives Identified	Concentration	LOQ*
GC/MS on derivatized ethanol extract	Ethanedioic acid, dibutyl ester (Oxalic)	0.16 ppm	0.01 ppm
	Propanedioic acid, dibutyl ester (Malonic)	0.10 ppm	0.01 ppm
	Dodecanoic acid, butyl ester (Lauric)	0.23 ppm	0.025 ppm
	Butanedioic acid, butyl ester (Succinic)	0.17 ppm	0.01 ppm
	Tetradecanoic acid, butyl ester (Myristic)	0.09 ppm	0.025 ppm
	Hexadecanoic acid, butyl ester (Palmitic)	0.61 ppm	0.01 ppm
	Octadecanoic acid, butyl ester (Stearic)	1.23 ppm	0.05 ppm
	1,4-Benzenedicarboxylic acid, bis(2-methylpropyl) ester	2.24 ppm	0.025 ppm
* Limit of quantification			
Test	Compounds Provisionally Identified	Molecular Weight	
LC/MS on ethanol extract	2-Ethylhexanoic acid	144.21	
	Lauric acid	200.32	
	Myristic acid	228.37	
	Palmitic acid	256.42	
	Stearic acid	284.48	

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
Decision tree for risk analysis of extractables and leachables

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Ding W. Martin J. Implementation of single-use technology in Biopharmaceutical Manufacturing. Bio Process Int 2009 (May): 46-51.