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Compound	Analytical Method	Abbreviation
Organic components	Total organic carbon analysis	TOC
Conducting ions	Conductivity analysis	σ
Acidic or basic compounds	Acidity or basicity measurement	pH
Acetate and formate	Ion-exchange liquid chromatography	IC
Nonvolatiles	Nonvolatile residue analysis	NVR
Functional groups of organic compounds	Fourier-transform infrared spectroscopy	FTIR
Chromophoric groups of organic compounds	Ultraviolet spectroscopy	UV
Volatile organics	Head-space gas chromatography with mass spectrometry	GC/MS
Semivolatiles	Direct-injection gas chromatography with mass spectrometry	GC/MS
Heat-sensitive organics; antioxidants	High-performance liquid chromatography with UV spectroscopy	HPLC/UV
Organic acids	Liquid chromatography with mass spectrometry	LC/MS
	Derivatization gas chromatography with mass spectrometry	GC/MS
Inorganic elements	Inductively coupled plasma mass spectrometry	ICP/MS

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
Tubing Properties		Extraction Conditions	
Material	SEBS thermoplastic elastomer*	Model solvents	DI water (18 M Ω cm); 100% ethanol (HPLC grade)
Liquid contact length	~4.5 m (15 ft)		
Weight	113 g/m (0.076 lb/ft)		
Wall thickness	3.2 mm (0.125 in)	Total volume	306 mL
Internal diameter	9.5 mm (0.75 in)	Time	72 hours
Liquid contact area	299 cm ² /m (14.1 in ² /ft)	Temperature	21 °C
Pretreatment	Gamma irradiation (50 kGy)	Agitation method	Orbital shaker

*Tubing base polymers: styrene-ethylene-butylene-styrene copolymer and polypropylene

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Sample	TOC	Conductivity	pH ¹	Acetate	Formate
Tubing extract	5.37 ppm	7.09 μ S/cm	5.6	0.034 ppm	0.032 ppm
Negative control	<LOD ²	1.06 μ S/cm	6.5	<LOD ²	<LOD ²

¹The pH measurement was performed separately using water for injection (WFI) at pH 6.6.

²Limit of detection (LOD): 0.002 ppm TOC; 0.003 ppm acetate; 0.001 ppm formate

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	Water Extract	Ethanol Extract
Volatiles ¹	2-methoxy-2-propanol	2,3,4-trimethylpentane
Semivolatiles ²	No peaks	Low MW aliphatic hydrocarbons; 1,3-ditertbutylbenzene (1,3-DTBB); 2,4-ditertbutylphenol (2,4-DTBP); 1-tridecanol; lauryl acrylate
Nonvolatiles and heat-sensitive compounds ³	No peaks	Irgafos* antioxidant and its degradants
Organic acids ⁴	Oxalic acid	Myristic, palmitic, and stearic acid
Inorganic elements ⁵	Na, Fe, Zn, and K <10 ppb; Ca 53 ppb	Na, Al, K, and Ca no higher than for the negative control

¹ LOD: 0.001–0.1 ppm; LOQ: 0.01–0.5 ppm

² Limit of detection: 0.01 ppm

³ Limit of detection: 0.1 – 1 ppm

⁴ LOD: 0.003–0.01 ppm; LOQ: 0.025 ppm

⁵ LOD: 0.01–0.14 ppb for water; 1–5 ppb for ethanol

* Irgafos is a proprietary stabilizer and trademark of Ciba Holding AG in Basel, Switzerland.


جدول ۶ - مشخصات و نتایج بررسی وضعیت ظروف			
Biocontainer Properties		Extraction Conditions	
Fluid contact layer material	Low-density polyethylene	Model solvents	DI water (18 MΩ·cm); 100% ethanol (HPLC grade)
Internal gas barrier layer material	Ethyl vinyl alcohol		
Outer mechanical barrier layer material	Low-density polyethylene		
Capacity	0.5 L	Total volume	2 L
Liquid contact area	516 cm ² (80 in ²)	Time	30 days
Pretreatment	Gamma irradiated 50 kGy	Temperature	40 °C

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¹ pH measurement was performed separately using WFI (water for injection) at pH 6.6

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	Water Extract	Ethanol Extract
Volatiles ¹	2-methyl-2-propanol; butanal; hexanal; 2-octanone	2-methylpentane; hexane; trimethylpentane; 3-methylheptane; 1-octene; octane
Semi-volatiles ²	No peaks	Low MW oligomers of polyethylene; 1,3-DTBB; 2,4-DTBP; 2-octanone; 1-heptadecanol; 1-octadecanol
Organic fatty acids ³	Oxalic acid 0.07 ppm	Succinic, palmitic, and stearic acid all below 0.06 ppm
Inorganic elements ⁴	Na, B, Mg, K, Ca, and Ba all at ppb levels	B 8.3 ppb

¹ LOD: 0.001-0.1 ppm; LOQ: 0.01-0.5 ppm² LOD: 0.01 ppm³ LOD: 0.003-0.01 ppm; LOQ: 0.025 ppm⁴ LOD: 0.01-0.14 ppb water; 1-5 ppb ethanol

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Ding W. Martin J. Implementation of single-use technology in Biopharmaceutical Manufacturing. Bioprocess Int 2009; 6: 45-50.